

Introducing NEW Microscopes

from  **ken-a-vision**
KNOWLEDGE THROUGH VISION

**FREY**
scientific

Take advantage of Frey Scientific's reduced prices on **NEW** microscopes from Ken-A-Vision.

Microscopes are a great tool for any classroom. They allow students to explore objects more in depth than with the human eye. Now with digital microscopes, teachers can help students see details easier than ever before.

You and your students can save the images and view them on a computer monitor. Students also can use the images to create more interesting presentations, lab reports or to demonstrate time progression in their findings.

Now at reduced prices!

Hurry. These special discount prices end July 1, 2010!



Introductory offer
25%
off list price



Digital Comprehensive Scope 2

Make it easy for students to share their investigations. Professional scope with sophisticated digital capabilities allows them to instantly display a specimen on a computer screen or through a projector, and to incorporate high-resolution images in documents or Web pages. PC and Mac compatible Applied Vision™ software captures single images, time-lapse sequences and movies, or digital movies at up to 30 fps. The instrument's modern, modular design features a stable base, built-in carry handle, and integrated accessory storage. USB cable included. Ten-year microscope warranty; one-year camera warranty.

Objectives: Achromatic 4X, 10X, 40XR, 100XR

Illumination: LED with rheostat

Stage: True mechanical

Focal Adjustments: Coaxial course and fine

Cameras: 1.3 MP CMOS (1280 x 1024)

Software: Applied Vision™ (PC and Mac)

Part #	Type	Reg. Price	Special Price
D1324387-628	Monocular	\$899.00	\$675.00
D1324388-628	Binocular	\$1299.00	\$975.00

D-Drop Ship

Save up to 25% on Comprehensive Ken-A-Vision Microscopes

Comprehensive Scope 2

A professional instrument priced to fit an educational budget, this scope sports a modern, modular design on a very stable base. Operates anywhere you need it on rechargeable NiMH batteries that last up to 40 hours on an 8-hour charge, or the AC adaptor, which conveniently stores in the unit. Features patented LED bright white light with long bulb life, reverse nosepiece for easy slide insertion and manipulation, and built-in carry handle. Ten-year warranty.

Objectives: Achromatic 4X, 10X, 40XR, 100XR

Eyepiece: 10x widefield with pointer

Illumination: LED with rheostat

Stage: Plain or True Mechanical (Bino only)

Focal Adjustments: Coaxial course and fine

Weight: 10.4 lbs

Height: 14"

Part #	Type	Reg. Price	Special Price
D1324384-628	Monocular	\$399.00	\$299.00
D1324385-628	Dual-Opposing	\$449.00	\$336.00
D1324386-628	Binocular	\$849.00	\$636.00

D-Drop Ship



D1324384



D1324386



Frey Scientific has created a "Guide to Digital Microscopy" to show teachers how easy it is to incorporate digital microscopy into their lessons. To learn more call 1-800-225-3739

Part #	Title	Reg. Price
N1327426-628	Guide to Digital Microscopy (Printed Binder)	\$99.95
N1327773-628	Microscopy CD-ROM (Individual)	\$99.95
N1329742-628	Microscopy CD-ROM (Network Version)	\$499.95

N-No Discounts Apply

A shipping and handling charge of 12% will be added to all orders, with a minimum of \$5.00

Contact your Frey Scientific Sales Representative at 1-800-225-3739
www.FreyScientific.SchoolSpecialty.com